

PRODUCT CODE CZ009

INTENDED USE


This reagent is intended for *in vitro* quantitative determination of Creatine Kinase (CK) in serum or plasma.

CLINICAL SIGNIFICANCE

Creatine Kinase is a cellular enzyme with wide tissue distribution in the body. Its physiological role is associated with adenosine triphosphate (ATP) generation for contractile or transport systems. Elevated CK values are observed in diseases of skeletal muscle and after myocardial infarction. Clinical diagnosis should not be made on a single test result; it should integrate clinical and other laboratory data.

PRINCIPLE

Creatine Kinase (CK) catalyses the reversible transfer of a phosphate group from phosphocreatine to ADP. This reaction is coupled to those catalyzed by hexokinase (HK) and glucose-6-phosphate dehydrogenase (G6P-DH).



+ REAGENT COMPOSITION

REAGENT 1

Imidazole,	125 mmol/L
pH 6.70	
D-Glucose	25 mmol/L
N-Acetyl-L-Cysteine	25 mmol/L
Magnesium acetate	12.5 mmol/L
NADP	2.4 mmol/L
EDTA	2.0 mmol/L
Hexokinase	6000 U/L

REAGENT 2

ADP	15.0 mmol/L
AMP	25 mmol/L
P1, P5-di-Adenosine-5-	102 µmol /L
pentaphosphate Glucose-6-phosphate	8000 U/L
dehydrogenase Creatinine phosphate	250 mmol/L

PREPARATION OF WORKING REAGENT

Mix 4 volumes of R1 with 1 volume of R2, Stability: 2 weeks at 2-8°C or 48 hours at room temperature (15-25°C)

STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C, protected from light and contaminations prevented during their use. **Signs of reagent deterioration:**

- Presence of particle and turbidity
- Blank absorbance (A) at 340nm ≥ 1

SPECIMEN

Fresh serum free of hemolysis or heparinized plasma

Stability 7 days at 2-8° C, protected from light. The Creatine Kinase activity decreases 10% after 1 day at 2-5° C or after 1 hour at 15-25°

C. PRUCEDURE

APPLICATION ON AUTO BIOCHEMISTRY ANALYSERS

PERFORMANCE CHARACTERISTICS

detection limit: 9.2 U/L = 153 nkat/L

linearity limit of 1300 U/L = 21671 nkat/L

for higher values dilute sample 1/2 with distilled water and repeat measurement.

– Repeatability (within

run): Mean Concentration

	C.V	n
175 U/L = 2917 nkat/	1.8 %	20
L 567 U/L = 9452 nkat/L	0.7 %	20

– Reproducibility (run to

run): Mean Concentration

	C.V	n
175 U/L = 2917 nkat/	1.3 %	25
L 567 U/L = 9452 nkat/L	1.1%	25

INTERFERENCES

Bilirubin (< 20 mg/dL) and hemoglobin (< 10 g/L) do not interfere. Lipemia (triglycerides > 5 g/L) interfere. Other drugs and substances may interfere



QUALITY CONTROL

It is recommended to use CK NAC control sera of known value.

If control values are found out of the defined range, check the instrument, reagents and technique for problems.

NORMAL VALUES

	U/L	nkat/L
Men, up to	38-174	633-2900
Women, up to	26-140	433-2334

Each laboratory should establish its own normal range representing patient population

SYMBOL ON LABELS

IVD in vitro diagnostics	 manufacturing date
LOT lot number	 expiry date
REF catalogue number	 manufacturer
 temperature limit	 instruction for use

BIBLIOGRAPHY

1. IFCC primary reference procedures for the measurement of catalytic activity concentrations of enzymes at 37°C, Part 2. Reference procedure for the measurement of catalytic concentration of creatine kinase. Clin Chem Lab Med 2002;40:635-642.
2. IFCC reference procedures for measurement of catalytic concentrations of enzymes: corrigendum, notes and useful advice. Clin Chem Lab Med 2010; 48: 615-621.
3. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics, 4th ed. Burtis CA, Ashwood ER, Bruns DE. WB Saunders Co, 2005.
4. Young DS. Effects of drugs on clinical laboratory tests, 5th ed. AACC Press, 2000.
5. Friedman and Young. Effects of disease on clinical laboratory tests, 4th ed. AACC Press, 2001..

