

PRODUCT CODE CS022

INTENDED USE

This reagent is intended for in vitro quantitative determination of total Bilirubin in serum or plasma.

CLINICAL SIGNIFICANCE

Bilirubin is a waste product derived from the heme moiety of the hemoglobin released from senescent or damaged erythrocytes that are destroyed in the reticuloendothelial cells. After production, bilirubin is transported to the liver in association with albumin. Inside the hepatocytes bilirubin is conjugated with glucuronic acid and it is excreted into bile. A number of inherited and acquired diseases affect production, uptake, metabolism, and excretion of bilirubin, resulting in hyperbilirubinemia^{1,2}. Unconjugated hyperbilirubinemia is seen in newborns (physiological jaundice), in increased red cell destruction (hemolytic anemia, extensive hematoma), in ineffective erythropoiesis and in some rare genetic diseases (Gilbert's syndrome, Crigler-Najjar syndrome). Conjugated hyperbilirubinemia is associated to a decreased excretion of bile due to liver diseases (hepatitis or cirrhosis) or to intrahepatic or extrahepatic cholestasis. Jaundice is a clinical manifestation of hyperbilirubinemia, consisting of deposition of bile pigments in the skin, resulting in a yellowish staining of the skin and mucous membranes. Clinical diagnosis should not be made on the findings of a single test result, but should integrate both clinical and laboratory data.

PRINCIPLE

Direct bilirubin in the sample reacts with 3,5-dichlorophenyl diazonium salt forming a colored complex that can be measured by spectrophotometry at 546 nm³. Both direct and indirect bilirubin couple with diazo in the presence of cetrimide^{4,5}. The terms "direct" and "total" refer to the reaction characteristics of serum bilirubin in the absence or presence of solubilizing (accelerating) reagents. The "direct" and "indirect" bilirubin are only approximately equivalent to the conjugated and unconjugated fractions.

CONTENTS:

Reagent 1: Hydrochloric acid 170 mmol/L, cetrimide 40 mmol/L, pH 0,9.

Reagent 2: 3,5-dichlorophenyl diazonium 1.5 mmol/L.

REAGENT PREPARATION

All reagents are ready to use.

REAGENT STORAGE AND STABILITY

Store at 2-8°C. Components are stable once opened until the expiry date marked in the label if they are stored at the recommended temperature, well closed and care is taken to prevent contamination during their use.

SPECIMEN

Serum and plasma collected by standard procedures. EDTA may be used as anticoagulants. Bilirubin concentration in serum and plasma is stable for 1 days at 20-25°C, 3 days at 4-8°C and 6 months at -20°C if protected from light

PRECAUTION

To avoid contamination, use clean laboratory wares. Avoid direct exposure of reagent to light.

APPLICATION ON AUTO BIOCHEMISTRY ANALYZERS

BR-160 AUTO BIOCHEMISTRY ANALYZER					
TEST NAME	BIL.T	SAMPLE VOLUME		5	
FULL NAME	BILIURUBIN TOTAL	REAGENT VOLUME		R1	R2
WAVE LENGTH	PRI. 546	SEC. 700	200	50	
ASSAY / POINT	1 POINT END	START 1	END 25		
DECIMAL PLACE	2				
UNIT	mg/dl				
LINEARITY RANGE	LOW 0.201	HIGH 38			

BR-120 AUTO BIOCHEMISTRY ANALYZER					
TEST NAME	BIL.T		SAMPLE VOLUME	5	
FULL NAME	BILIURUBIN TOTAL				
	PRI.	SEC.	REAGENT VOLUME	R1	R2
WAVE LENGTH	546	700		200	50
ASSAY /POINT	1 POINT END	START	END		
		1	39		
DECIMAL PLACE	2				
UNIT	mg/dl				
LINEARITY RANGE	LOW	HIGH			
	0.201	38			

- Detection limit: 0.201 mg/dL = 3.43 μ mol/L.
- Linearity: 38 mg/dL = 650 μ mol/L.
- Precision:

Mean concentration	Repeatability (CV)	Within-laboratory (CV)
2.10 mg/dL = 36.0 μ mol/L	2.9 %	6.1 %
4.72 mg/dL = 80.6 μ mol/L	1.9 %	4.7 %

Trueness: Results obtained with this reagent did not show systematic differences when compared with reference reagents. Details of the comparison experiments are available on request.

NORMAL RANGE

Adult: Up to 2.0 mg/dL = 34 μ mol/L

L Newborns:

Age	premature	full-term
Up to 24 h	1.0-8.0 mg/dL = 17-137 μ mol/L	2.0-6.0 mg/dL = 34-103 μ mol/L
Up to 48 h	6.0-12.0 mg/dL = 103-205 μ mol/L	6.0-10 mg/dL = 103-171 μ mol/L
3-5 days	10-14 mg/dL = 171-239 μ mol/L	4.0-8.0 mg/dL = 68-137 μ mol/L

QUALITY CONTROL

It is recommended to use the Biochemistry Control Serum Biotrol level I & Biotrol level

II

LIMITATION & PRECAUTIONS:

– Interferences: hemolysis (hemoglobin 500 mg/dL) do not interfere. Lipemia (triglycerides 1300 mg/dL) interfere. Other drugs and substances may interfere.

SYMBOL ON LABELS

IVD in vitro diagnostics	manufacturing date
LOT lot number	expiry date
REF catalogue number	manufacturer
temperature limit	instruction for use

BIBLIOGRAPHY

1. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics, 4th ed. Burtis CA, Ashwood ER, Bruns DE. WB Saunders Co, 2005.
2. Friedman and Young. Effects of disease on clinical laboratory tests, 4th ed. AACC Press, 2001.
3. Thaler M, Luppa PB and Schlebusch H. Bilirubin measurement – an updated survey. J Lab Med 2008; 32:1-9.
4. Zoppi F, Peracino A, Fenili D, Marcovina S and Ramella C. Metodo per la determinazione della bilirubina totale e coniugata. Uso di un tensioattivo cationico come agente solubilizzante. Giorn It Chim Cl 1976; 1:343-359.
5. Pearlman FC and Lee RTY. Detection and measurement of total bilirubin in serum, with use of surfactants as Solubilizing agents. Clin Chem 1974; 20: 447-453.
6. Young DS. Effects of drugs on clinical laboratory tests, 5th ed. AACC Press, 2000.

