

PRODUCT CODE CS027

INTENDED USE



This reagent is intended for *in vitro* quantitative determination of Urea in serum, plasma or urine.

CLINICAL SIGNIFICANCE

Urea is the end product of protein metabolism. It is synthesized in the liver from the ammonia produced by the catabolism of amino acids. It is transported by the blood to kidneys from where it is excreted. Increased level associated in renal disease, urinary obstruction, shock, congestive heart failure and burns. Decreased levels are found in liver failure and pregnancy.

PRINCIPLE

Urea in the sample hydrolyzed enzymatically into ammonia (NH₄) and carbon dioxide (CO₂). The ammonia formed further combined with α-ketoglutarate in a reaction catalyzed by glutamate dehydrogenase (GLDH) with simultaneous oxidation of NADH to NAD⁺



The decrease in concentration of NADH, is proportional to urea concentration in the sample¹.

REAGENT

COMPOSITION

REAGENT 1, Buffer	80 mmol/L
α-Ketoglutarate	6 mmol/L
Urease	75000 U/L
REAGENT 2	
(Enzymes) GLDH	60000 U/L
NADH	0.32 mmol/L

REAGENT PREPARATION

R1 and R2, ready-to-use and stable upto the expiry date if contamination is avoided and stored at 2-8°C and protect from light. **Signs of reagent deterioration**

Presence of turbidity

Blank absorbance (A) at 340 nm < 1.00

SPECIMEN

Serum, heparinized plasma; don't use ammonium salts or fluoride as anti

coagulant. Urine sample 1+49 with distilled water before the assay, multiply the results by 50. Preserve urine sample at pH<4.

Urea is stable for 5 days when stored at 2-8°C.

PROCEDURE

APPLICATION ON AUTO BIOCHEMISTRY ANALYSERS

BR-160 AUTO BIOCHEMISTRY ANALYZER					
TEST NAME	UREA		SAMPLE VOLUME	3	
FULL NAME	UREA uv				
	PRI	SEC	REAGENT VOLUME	R1	R2
WAVE LENGTH	340	NONE		240	60
ASSAY /POINT	fixed-time	START	END		
		14	20		
DECIMAL PLACE	2				
UNIT	mg/dl				
LINEARITY RANGE	LOW	HIGH			
	1	400			

BR-120 AUTO BIOCHEMISTRY ANALYZER					
TEST NAME	UREA		SAMPLE VOLUME	3	
FULL NAME	UREA ur.				
	PRI.	SEC.	REAGENT VOLUME	R1	R2
WAVE LENGTH	340	NONE		240	60
ASSAY POINT	fixed-time	START	END		
		22	31		
DECIMAL PLACE	2				
UNIT	mg/dl				
LINEARITY RANGE	LOW	HIGH			
	2	400			

QUALITY CONTROL

Control Sera are recommended to monitor the performance of assay procedure. If control values are found out of the range, check the instrument, reagent, and standard for problems. Each laboratory should establish its own Quality Control scheme and corrective actions if the controls don't meet the acceptable tolerances.

REFERENCE VALUES ⁴

Serum or plasma: 15 – 45 mg/dL = 2.5 – 7.5 mmol/L

Urine: 26 – 43g/24 h = 428 – 714 mmol/24 h

These values are for orientation purpose; each laboratory should establish its own reference range.

PERFORMANCE CHARACTERISTICS

Measuring Range: From detection limit 1 mg/dL to linearity limit 400 mg/dL.

If the concentration is greater than linearity limit dilutes 1: 2 the sample with normal saline and multiply the result by

2. Precision

Intra-assay n=20	Mean (mg/dL)	SD	% CV
Level I	40.70	0.88	2.16
Level II	130	1.02	0.78

Inter-assay n=20	Mean (mg/dL)	SD	% CV
Level I	40.50	1.19	2.94
Level II	128	2.07	1.61

Sensitivity/Limit of Detection

The lower limit of detection is 1 mg/dL.

Accuracy

Results obtained using Bio Research reagent (y) didn't show any systematic differences when compared with other commercial reagent (x). The results obtained using 50 samples was the following,

Correlation coefficient: (r)² :0.998

Regression equation y= 1.5759x – 1.1577

The results of performance characteristics depend on the analyzer used.

INTERFERENCES

It is recommended to use heparin as anticoagulant. Do not use ammonium salts or fluorides¹.

A list of drugs and other interfering substances with urea determination has been reported by Young et al².

3. NOTES

- Glassware and distilled water must be free from ammonia and ammonium salt.
- Calibration with aqueous standard may cause a systematic error in automatic procedure. In these cases, it is recommended to use a serum calibrator.
- Use clean disposable pipette tips for dispensing.
- Bio Research has instruction sheets for several automated analyzers. Instructions for many of them are available on request.

SYMBOL ON LABELS

IVD	in vitro diagnostics		manufacturing date
LOT	lot number		expiry date
REF	catalogue number		manufacturer
	temperature limit		instruction for use

BIBLIOGRAPHY

- Kaplan A. Urea. Kaplan A et al. Clin Chem the C.V. Mosby Co. St Louis, Toronto, Princeton 1984; 1257-1260 and 437 and 418.
- Young DS. Effects of drugs on Clinical Lab. Tests, 4th ed AACC Press, 1995.
- Young DS. Effects of disease on Clinical Lab. Tests, 4th ed AACC 2001.
- Burtis A et al. Tietz Textbook of Clinical Chemistry, 3rd ed AACC 1999.
- Tietz N W et al. Clinical Guide to Laboratory Tests, 3rd ed AACC 1995.