

PRODUCT CODE CZ005

INTENDED USE

This reagent is intended for *in vitro* quantitative determination of AST/GOT in serum or plasma.

METHOD


Optimized UV KINETIC-TEST according to IFCC (International Federation of Clinical chemistry)

CLINICAL SIGNIFICANCE

The Aspartate aminotransferase (AST/GOT) a cellular enzyme, it is present in most of the tissues. Especially in cardiac muscle, liver cells, skeletal muscle & kidneys. Injury to these tissues results in the release of the enzyme in blood stream. Increased levels are found in myocardial infarction. The duration & extent of increase is related to the infarct. GOT determination is of considerable value to differentiate myocardial infarction from other cardiac disorders. Increased levels are also found in various types of liver disease, skeletal muscle trauma & in renal diseases. Decreased levels may be found in pregnancy, Beriberi & Diabetic ketoacidosis.

PRINCIPLE

the (AST/GOT) catalyzes the transfer of the amino group from L-aspartate to 2-Oxoglutarate to yield oxaloacetate and L-glutamate. The oxaloacetate undergoes reduction with simultaneous oxidation of NADH to NAD⁺ in the malate dehydrogenase (MDH) catalyzed indicator reaction. The resulting rate of decrease in absorbance at 340nm is directly proportional to the AST activity. Lactate dehydrogenase (LDH) is added to prevent interference from endogenous pyruvate which is normally present in serum.



REAGENT COMPOSITION

REAGENT 1 (ENZYME)

REAGENT) Tris pH 7.8	80 mmol/L
L-Aspartate	240 mmol/L
MDH	> 600 U/L
LDH	≥ 600 U/L

REAGENT 2

(SUBSTRATE) 2-

Oxoglutarate	12 mmol/L
PNADH	0.18 mmol/L

REAGENT

PREPARATION Reagent is ready for use.

SPECIMEN

Serum, heparinized plasma

PRUCEDURE

APPLICATION ON AUTO BIOCHEMISTRY ANALYSERS

BR-160 AUTO BIOCHEMISTRY ANALYZER					
TEST NAME	AST		SAMPLE VOLUME	20	
FULL NAME	Aspartate Amino Transferase				
	PRI	SEC	REAGENT VOLUME	R1	R2
WAVE LENGTH	340	NONE	200	50	
ASSAY /POINT	KINETIC	START	END		
		14	22		
DECIMAL PLACE	2				
UNIT	u/L				
LINEARITY RANGE	LOW	HIGH			
	0	467			

BR-120 AUTO BIOCHEMISTRY ANALYZER					
TEST NAME	lipase		SAMPLE VOLUME	4	
FULL NAME	lipase				
	PRL	SEC.	REAGENT VOLUME	R1	R2
WAVE LENGTH	578	700		200	50
ASSAY /POINT	FIXED - TIME	START	END		
		26	34		
DECIMAL PLACE	2				
UNIT	u/L				
LINEARITY RANGE	LOW	HIGH			
	0	300			

PRECAUTION

- The reagents contain sodium azide as preservative. Do not swallow and avoid contact with skin and mucous membranes.
- To avoid contamination, use clean laboratory wares. Avoid direct exposure of reagent to light.

Measuring range:

Up to linearity limit of 467 U/L.

If the results obtained were greater than linearity limit, dilute the sample 1/10 with NaCl 9 g/L and multiply the result by 10.

NORMAL RANGE

Men up to 37 U/L

Women up to 31 U/L

Each laboratory should establish reference ranges for its own patients' population.

QUALITY CONTROL

All control sera with values determined by this method can be used.

SYMBOL ON LABELS

IVD in vitro diagnostics	 manufacturing date
LOT lot number	 expiry date
REF catalogue number	 manufacturer
 temperature limit	 instruction for use

BIBLIOGRAPHY

- Clin. Chem. ACTA 105 (1980) S. 147-172 Synopsis Der Leberkrankheiten : H. Wallhofer, E. Schmidt.
- .Thefeld W. ET. AI. DT . MED. WSCHR. 99 (1974) 343.

